- 18- A hysteroscopic tubal sterilization using flexible micro inserts (Essure) is performed, e procedure time is about 30 minutes (from insertion to removal of hysteroscope), and there was a concern regarding possible perforation due to a sudden loss of resistance at insertion. What is the suggested method of confirming that the procedure has been successful?

  Options
  - A. Hysterosalpingogram.
  - B. Magnetic resonance imaging.
  - C. Transvaginal ultrasound scan.
  - D. Pelvic x-ray.
  - E. Diagnostic laparoscopy.
- 19- An 18-year-old nulliparous patient attends the gynaecology clinic for heavy painful menstrual bleeding. She is in a sexual relationship with a new partner and both have tested negative for STIs, though she has had a previously treated chlamydia infection. She does not want any children in the near future. She takes lamotrigine for epilepsy control and has not suffered any seizures for 3 years. General and pelvic examinations are normal, and her BMI is 30. What is the best appropriate method of contraception in her situation? Options
  - A. Combined oral contraceptive pill.
  - B. Depot medroxyprogesterone acetate.
  - C. Evra patch.
  - D. Levonorgestrel-releasing intrauterine system.
  - E. NuvaRing
- 20- A 35 year-old nulliparous patient is HIV positive and takes efavirenz and nevirapine. She is using a 30-mcg estrogen combined oral contraceptive preparation as well as barrier contraception. She finished a pill packet 10 days ago, forgot to restart again, and had unprotected sexual intercourse 2 days ago. She also missed two pills in the first week of the previous pill packet. A chlamydia urine polymerase chain reaction test is positive. She attends asking for emergency contraception. What is the most appropriate method of emergency contraception in this situation? Options
  - A. 1.5 mg levonorgestrel as soon as possible.
  - B. 3 mg levonorgestrel as soon as possible.
  - C. LNG-IUS as soon as possible.
  - D. e copper IUD as soon as possible.
  - E. Ulipristal acetate single 30 mg tablet as soon as possible

